



An Innovative Sampling and Solventless Field Extraction Technique for Determination of Explosives Residues in Groundwater

**David E. Splichal, Denise K. MacMillan
and P. Stephen Schnitker
Environmental Laboratory, Omaha, NE**

John P. Shannon, ASI, Vicksburg, MS

Bradley Varhol, EON Products, Inc., Snellville, GA


**46th Rocky Mountain Conference on Analytical Chemistry
August 3, 2004**

Long Term Monitoring Focus Area

- Long Term Monitoring (LTM) of groundwater :
 - Required component of closure.
 - All military services, other Federal agencies (e.g., DOE), states, and responsible parties share similar responsibility.
 - Costs associated with sampling and laboratory analysis over 10 years estimated to approach \$500M.
 - Sample collection and laboratory analysis
 - ◆ 70% of the total monitoring cost.
 - ◆ 50% of the total investigation cost.

Interim Improvements

- ✓ COT/GOTS
- ✓ QA Processes & Protocols
- ✓ Direct Push Wells & Samplers
- ✓ Solventless Extraction Technologies
- ✓ Solventless Extraction Technologies Interfaced to Miniature GC






Solventless Extraction Technologies

US Army Engineer Research & Development Center



- ✓ Identify & Develop Solventless Ext Technologies
- ✓ Perform Lab & Field Studies
- ✓ Investigating use of Twister and SPME for ORCs
- ✓ Evaluation of Potential for On-Site Extraction



Fixed Lab – Sample Receipt

- Lots of coolers.
- Lots of ice.
- Lots of heavy lifting.
- Lots of sample bottles.
- Lots of time.
- Lots of MONEY! Approximately \$500 to ship 20 samples in seven coolers from field to the laboratory via overnight delivery services.



Fixed Lab – Extraction Procedures

- 2 x 1-L amber bottles/sample (44 bottles for 20 samples including QC samples)
- Double salting out procedure – SW-846, method 8330, 1994-1997. Preparation time: 6-8 hours
- Solid Phase Extraction – SW-846, method 3535A, Jan. 1998. Preparation time: 3-5 hours.





Solventless Extraction

US Army Engineer Research & Development Center

- **Collect known volume of sample.**
- **Perform Solid Phase Extraction in the field.**
- **Send only SPE cartridges to the laboratory for analysis.**

**ULTIMATE GOAL: SAVE MONEY THROUGH
LOWER SHIPPING COSTS.**



InSTED

ECB prototype

“*In-Situ* Tubular Extraction Device”

- Uses a modified HydraSleeve™ for sample collection (EON Products, Inc, Snellville, GA).
- SPE cartridge (Porapak RDx, Waters Corporation, Milford, MA).
- Submersible pump assembly (Clark Solutions, Hudson, MA).




InSTED

US Army Engineer Research & Development Center



ECB prototype *In-Situ* Tubular Extraction Device



InSTED

Experimental Design

- Blank Spikes Study
- Real World Compositated Samples Study
- Initial Time Study
- Standpipe Study



InSTED

Results for Blank Spike Extraction of Explosives at 10 ug/L (ppb)

InSTED

Traditional Method 8330

Analyte

% Recovery

Recovery

HMX

82

90

RDX

93

99

Tetryl

106

89

TNT

93

89

2,4-DNT

88

83

US Army Engineer Research & Development Center





InSTED

Results for Composited Real World Extraction of Explosives

<u>Analyte</u>	<u>InSTED</u> <u>ug/L (ppb)</u>	<u>Traditional Method 8330</u> <u>ug/L (ppb)</u>
HMX	21.4	21.6
RDX	153	158
4-A-DNT	1.5	1.5
2-A-DNT	0.1	0.1

US Army Engineer Research & Development Center





InSTED

Results for Composited Real World Extraction of Explosives – Time Study

US Army Engineer Research & Development Center

<u>Analyte</u>	InSTED	InSTED	<u>% Loss</u>
	immediate	10-days later	
	<u>ug/L (ppb)</u>	<u>ug/L (ppb)</u>	
HMX	21.4	19.3	10
RDX	153	137	10
4-A-DNT	1.5	1.0	33
2-A-DNT	0.1	0.06	40





InSTED

Results for Standpipe Extraction of Explosives at 10 ug/L (ppb)

<u>Analyte</u>	InSTED	Traditional Method 8330
	<u>% Recovery</u>	<u>Acceptable Recovery Ranges</u>
HMX	81	39 - 126
RDX	86	35 - 119
Tetryl	43	14 - 120
TNT	85	71 - 117
2,4-DNT	89	76 - 110



Potential Cost Savings

- **Traditional Sampling Costs:**

\$120 sample bottles, \$560 seven coolers, \$150 technician time.

TOTAL: \$830.

- **InSTED Costs:**

\$400 for HydraSleeves, \$80 for one cooler.

TOTAL: \$480.

COST SAVINGS: 85% in shipping charges.

40% in total sampling costs.

Future Work

- **Complete Time Studies:** How long can the SPE cartridges be held on-site? Ice necessary?
- **Surrogate Spike:** Monitor extraction efficiency. Add to bag prior to sending to field?
- **Field Evaluation:** Side-by-side sampling at existing site(s).
- **Independent Laboratory Evaluation:** Send SPE cartridges to referee laboratory to compare results.
- **GOAL:** Incorporate into COE EM-200-1-3 or into SW-846.
- **GOAL:** Investigate patentability.

Conclusions

- Potential for use in field applications.
- Results from the InSTED compare favorably with results from fixed laboratory extraction.
- Cost savings are significant.

Acknowledgements

- **John M. Cullinane, , USACE – ERDC, EQT Program Manager**
- **William Houser, AEC, Focus Area Manager**
- **Douglas B. Taggart, USACE – ECB, Branch Chief**

- **Project: Solventless Extraction Technologies for Long Term Monitoring of Military Unique Compounds. Work Item Code: B730JF.**



This presentation represents the personal views of the authors and does not suggest or reflect the policy, practices, program or doctrine of the US Army or the Government of the United States. The contents are not to be used for advertising or promotional purposes. Citation of brand names does not constitute official endorsement or approval of the products.